

Technical specification

Configura® Comfort 2

	Small	Medium
Overall chair width	740 mm (29")	790 mm (31")
Overall height	1060 mm (41.7")	1110 mm (43.7")
Overall length	885 mm (34.8")	930 mm (36")
Overall chair depth with legs elevated	1565 mm (61.6")	1660 mm (65.4")
Internal seat width	350 mm - 500 mm (14" - 20")	400 mm - 550 mm (16" - 22")
Floor to seat height	400 mm (16")	450 mm (18")
Internal seat depth Position 1	450 mm (18")	450 mm (18")
Internal seat depth Position 2	500 mm (20")	500 mm (20")
Backrest height from seat	710 mm (28")	710 mm (28")
Chair to wall/furniture	570 mm (22.4")	630 mm (24.8")
Backrest angle - std	100°	100°
Backrest angle - recline	128°	128°
Rise angle - std	23°	23°
Rise angle - med rise	14°	14°
Rise angle - flat	2°	2°
Rise height at front of seat - std	510 mm (20.1")	545 mm (21.5")
Rise height at front of seat - med	570 mm (22.4")	620 mm (24.4")
Rise height at front of seat - flat	660 mm (26")	710 mm (28")
TIS angle	29°	29°
Standard seat angle in neutral	6°	6°
Forward tilt angle	n/a	n/a
Armrest height	160 mm (6.3")	180 mm (7.1")
Weight capacity	160 kg (25 st)	160 kg (25 st)
Overall weight	65 kg (143 lbs)	68 kg (150 lbs)