

# Ascenda<sup>®</sup>

Transfer from lying to standing.  
Effortlessly.

The **Ascenda** is a 'riser-bed' designed to facilitate effortless transfers from lying to standing, enabling individuals to maintain dignity and mobility at home while reducing care package requirements.



**Tilt-in-space mechanism**



**Sleeping position**



**Rise function**



**Depth-adjustable**



**Pillow clips**



**Handset**

Overall dimensions	2150 mm (L) x 900 mm (W)
Mattress size	2100 mm (L) x 900 mm (W) x 150 mm (D)
Safe working load	245 kg
Patient weight	225 kg
Liquid ingress protection	IPX6
Overall weight	75 kg
Part number	CHB-0-CHI-000
25 st mattress part number	VISMAT25-CHI-FM5-000
35 st mattress part number	VISMAT35-CHI-FM5-000