


FloorBed¹ Plus

A wider floor-level nursing bed


Designed for the plus size user, FloorBed 1 Plus is 15 cm wider than a standard 90 cm wide bed, giving the plus size user more space without going to the full width of a bariatric bed.



Additional features:



Ultra Low function



Raises to 650 mm



Trendelenburg



Anti-Trendelenburg

Overall dimensions	1075 x 2255 mm (42.3" x 88.8")
Mattress size	1050 x 2000 mm (41.3" x 78.7")
Safe working load	260 kg (41 st, 573 lbs)
Patient weight	225 kg (35 st, 496 lbs)
Lifting range	78 - 650 mm (3.1" - 25.6")
Profiling functions	Kneebreak, Anti-Trendelenburg, Trendelenburg
Part number	NSB-0-FL5-100