


# Allevia<sup>®</sup>Harmony

This static mattress with in-built ventilation channels provides optimal pressure redistribution, without compromising on comfort.



|                        |                                             |
|------------------------|---------------------------------------------|
| <b>Sizes</b>           | 2000 x 900 x 150 mm<br>2000 x 1050 x 150 mm |
| <b>Weight capacity</b> | 160 kg (25 st)                              |
| <b>Warranty</b>        | 2 years                                     |
| <b>Cover</b>           | Four-way stretch PU                         |
| <b>Fire retardancy</b> | BS7177 Medium Hazard                        |
| <b>Seam type</b>       | Welded                                      |
| <b>Layered foams</b>   | ✓                                           |
| <b>Zip</b>             | ✓                                           |
| <b>Water resistant</b> | ✓                                           |