



Allevia® Comfort

A high quality mattress suitable for low to medium risk residents

This mattress has a base layer of high-density foam for support, with a top layer of soft, pressure-reducing visco foam. The cover is made of four-way stretch, vapour permeable fabric, for additional comfort.



Size	2000 x 900 x 150 mm (79 x 35 x 6") 2000 x 1050 x 150 mm (79 x 41 x 6")
Weight capacity	160 kg (25 st)
Warranty	2 years
Cover	Four-way stretch PU
Fire retardancy	BS EN 7177 Medium hazard (Ignition sources 0, 1 and 5)
Vapour permeable	✓
Layered foams	✓
Zip	✓
Water resistant	✓

