


Allevia® Comfort FirmEdge

Designed for additional support when side rails are used


The **FirmEdge** on the sides and ends of the mattress offers more support for transfers, and helps service users orient themselves in the bed. The **FirmEdge** also reduces the gaps between the mattress and head and foot board, to help reduce any risk of entrapment.



FirmEdge helps users to **orient themselves** in bed



Size	1950 x 900 x 150 mm (77 x 35 x 6")
Weight capacity	160 kg (25 st)
Warranty	2 years
Cover	Four-way stretch PU
Fire retardancy	BS EN 7177 Medium hazard (Ignition sources 0, 1 and 5)
Vapour permeable	✓
Layered foams	✓
Zip	✓
Water resistant	✓



Four-way stretch **vapour permeable** cover