




# Allevia® Prima

Featuring a castellated foam surface with a specialist heel management area

The specialist design and use of multi-stretch material enables the very best pressure re-distribution for most patient groups. The pressure-reducing castellated central zone maximises re-distribution in the most vulnerable areas, so that the mainstream patient weight can be managed effectively.



<b>Size</b>	2000 x 900 x 150 mm (79 x 35 x 6") 2000 x 1050 x 150 mm (79 x 41 x 6")
<b>Weight capacity</b>	168 kg (26 st)
<b>Warranty</b>	2 years
<b>Cover</b>	Four-way stretch PU
<b>Fire retardancy</b>	BS EN 7177 Medium hazard (Ignition sources 0, 1 and 5)
<b>Vapour permeable</b>	✓
<b>Layered foams</b>	✓
<b>Zip</b>	✓
<b>Water resistant</b>	✓



Four-way stretch **vapour permeable** cover